

## Practical Guide International EDEXCEL

This guide includes details about the core practicals for A-level chemistry. It also contains information about other experiments that often occur in A-level examinations. You may be asked to describe these experiments in details or be asked about reasons for doing individual steps.

You may be asked about other unfamiliar experiments but these will be using the skills and techniques that are described in the following experiments.

### Safety and hazards

- Irritant - dilute acid and alkalis- wear goggles
- Corrosive- stronger acids and alkalis wear goggles
- Flammable – keep away from naked flames
- Toxic – wear gloves- avoid skin contact- wash hands after use
- Oxidising- Keep away from flammable / easily oxidised materials

Hazardous substances in low concentrations or amounts will not pose the same risks as the pure substance.

### Measuring gas volumes


#### Using a gas syringe

Gas syringes can be used for a variety of experiments where the volume of a gas is measured, possibly to work out moles of gas or to follow reaction rates.

Potential errors in using a gas syringe


- gas escapes before bung inserted
- syringe sticks
- some gases like carbon dioxide or sulphur dioxide are soluble in water so the true amount of gas is not measured.

Make sure you don't leave gaps in your diagram where gas could escape



The volume of a gas depends on pressure and temperature so when recording volume it is important to note down the temperature and pressure of the room.

Alternatively gas volumes can be measured 'over water' with an up-turned measuring cylinder in a trough of water




Moles of gas can be calculated from gas volume (and temperature and pressure) using ideal gas equation  $PV = nRT$  or using the molar gas volume ( $1\text{mol gas} = 24\text{dm}^3$  at room temperature and pressure)

## Core Practical 1: Measure the molar volume of a gas

### Detailed method

1. Measure 30 cm<sup>3</sup> of 1 mol dm<sup>-3</sup> ethanoic acid and transfer to a conical flask.
2. Attach conical flask to gas syringe or use collection over water method (see previous page)
3. Measure the mass of a weighing bottle with approximately 0.05 g of calcium carbonate
4. Add the calcium carbonate to the conical flask- quickly resealing the bung so no gas escapes
5. Measure the final total volume of gas
6. Reweigh the empty weighing bottle test tube from step 3
7. Repeat the experiment several more times, increasing the mass of calcium carbonate by about 0.05 g each time.



### Analysis

From the graph read the volume of CO<sub>2</sub> given off with 0.25 g CaCO<sub>3</sub>

Work out the moles of CaCO<sub>3</sub> in 0.25g = 0.25/100.1 = 2.5 x 10<sup>-3</sup>

Assume the moles of CO<sub>2</sub> = moles of CaCO<sub>3</sub>

Work out molar volume of CO<sub>2</sub> = volume of CO<sub>2</sub>/ moles of CO<sub>2</sub>

### Method for using a gas syringe to calculate the Mr of propanone

1. Extract 0.20 cm<sup>3</sup> of propanone into a hypodermic syringe and then measure the mass of this syringe
2. using hand protection, remove a gas syringe from the oven and note the volume of air already in the barrel – about 5 cm<sup>3</sup>.
3. inject the propanone through the self-seal cap into the barrel. The plunger will move straight away.
4. Put the gas syringe back into the oven.
5. Measure the mass of the empty hypodermic syringe immediately.
6. After a few minutes measure the volume of the gas in the gas syringe, record the temperature of the oven shelf and the pressure of the room.

**Example 1** : 0.150g of a volatile liquid was injected into a sealed gas syringe. The gas syringe was placed in an oven at 70°C at a pressure of 100kPa and a volume of 80cm<sup>3</sup> was measured. Calculate the Mr of the volatile liquid (R = 8.31)

$$\text{moles} = PV/RT$$

$$= 100\,000 \times 0.00008 / (8.31 \times 343)$$

$$= 0.00281 \text{ mol}$$

$$Mr = \text{mass}/\text{moles}$$

$$= 0.15 / 0.00281$$

$$= 53.4 \text{ g mol}^{-1}$$


$$100 \text{ kPa} = 100\,000 \text{ Pa}$$

$$80 \text{ cm}^3 = 0.00008 \text{ m}^3$$

## Heating in a crucible

This method could be used for measuring mass loss in various thermal decomposition reactions and also for mass gain when reacting magnesium in oxygen.

The water of crystallisation in calcium sulfate crystals can be removed as water vapour by heating as shown in the following equation.




### Method.

- Weigh an empty clean dry crucible and lid .
- Add 2g of hydrated calcium sulfate to the crucible and weigh again
- Heat strongly with a Bunsen for a couple of minutes
- Allow to cool
- Weigh the crucible and contents again
- Heat crucible again and reweigh until you reach a constant mass ( do this to ensure reaction is complete).

Large amounts of hydrated calcium sulfate, such as 50g, should not be used in this experiment as the decomposition is likely to be incomplete.

The crucible needs to be dry otherwise a wet crucible would give an inaccurate result. It would cause mass loss to be too large as the water would be lost when heating.

The lid improves the accuracy of the experiment as it prevents loss of solid from the crucible but should be loose fitting to allow gas to escape.



Small amounts of the solid , such as 0.100 g, should **not** be used in this experiment as the percentage uncertainties in weighing will be too high.

**Example 2.** 3.51 g of hydrated zinc sulfate were heated and 1.97 g of anhydrous zinc sulfate were obtained. Use these data to calculate the value of the integer  $x$  in  $\text{ZnSO}_4 \cdot x\text{H}_2\text{O}$

Calculate the mass of  $\text{H}_2\text{O} = 3.51 - 1.97 = 1.54\text{g}$

Calculate moles of $\text{ZnSO}_4$	$= \frac{1.97}{161.5}$	Calculate moles of $\text{H}_2\text{O}$	$= \frac{1.54}{18}$
	$= 0.0122$		$= 0.085$

Calculate ratio of mole of $\text{ZnSO}_4$ to $\text{H}_2\text{O}$	$= \frac{0.0122}{0.0122}$		$= \frac{0.085}{0.0122}$
	$= 1$		$= 7$
			$X = 7$

## CORE PRACTICAL 2: To determine the enthalpy change of a reaction using Hess's Law

### Measuring the Enthalpy Change for a Reaction Experimentally

#### Calorimetric method

For a reaction in solution we use the following equation

**energy change = mass of solution x heat capacity x temperature change**

$$Q \text{ (J)} = m \text{ (g)} \times c_p \text{ (J g}^{-1}\text{K}^{-1}) \times \Delta T \text{ (K)}$$

This equation will only give the energy for the actual quantities used. Normally this value is converted into the energy change per mole of one of the reactants. (The enthalpy change of reaction,  $\Delta H$ )


#### Calorimetric method

One type of experiment is one in which substances are mixed in an insulated container and the temperature rise measured.

This could be a solid dissolving or reacting in a solution or it could be two solutions reacting together

##### General method


- washes the equipment (cup and pipettes etc) with the solutions to be used
- dry the cup after washing
- put polystyrene cup in a beaker **for insulation and support**
- Measure out desired volumes of solutions with volumetric pipettes and transfer to insulated cup
- clamp thermometer into place making sure the thermometer bulb is immersed in solution
- measure the initial temperatures of the solution or both solutions if 2 are used. Do this every minute for 2-3 minutes
- At minute 3 transfer second reagent to cup. If a solid reagent is used then add the solution to the cup first and then add the solid weighed out on a balance.
- If using a solid reagent then use 'before and after' weighing method
- stirs mixture (**ensures that all of the solution is at the same temperature**)
- Record temperature every minute after addition for several minutes



If the reaction is slow then the exact temperature rise can be difficult to obtain as cooling occurs simultaneously with the reaction

To counteract this we take readings at regular time intervals and extrapolate the temperature curve/line back to the time the reactants were added together.

We also take the temperature of the reactants for a few minutes before they are added together to get a better average temperature. If the two reactants are solutions then the temperature of both solutions need to be measured before addition and an average temperature is used.



##### Errors in this method

- energy transfer from surroundings (usually loss)
- approximation in specific heat capacity of solution. The method assumes all solutions have the heat capacity of water.
- neglecting the specific heat capacity of the calorimeter- we ignore any energy absorbed by the apparatus.
- reaction or dissolving may be incomplete or slow.
- Density of solution is taken to be the same as water.

Read question carefully. It may be necessary to describe:

- Method
- Drawing of graph with extrapolation
- Description of the calculation

## Calculating the enthalpy change of reaction, $\Delta H$ , from experimental data

### General method

1. Using  $q = m \times c_p \times \Delta T$  calculate energy change for quantities used
2. Work out the moles of the reactants used
3. Divide  $q$  by the number of moles of the reactant not in excess to give  $\Delta H$
4. Add a sign and unit (divide by a thousand to convert  $\text{J mol}^{-1}$  to  $\text{kJ mol}^{-1}$ )

The heat capacity of water is  $4.18 \text{ J g}^{-1}\text{K}^{-1}$ . In any reaction where the reactants are dissolved in water we assume that the heat capacity is the same as pure water.

Also assume that the solutions have the density of water, which is  $1 \text{ g cm}^{-3}$ . Eg  $25 \text{ cm}^3$  will weigh  $25 \text{ g}$

**Example 3.** Calculate the enthalpy change of reaction for the reaction where  $25.0 \text{ cm}^3$  of  $0.20 \text{ M}$  copper sulfate was reacted with  $0.01 \text{ mol}$  (excess of zinc). The temperature increased  $7.0^\circ\text{C}$ .

Step 1: Calculate the energy change for the amount of reactants in the calorimeter.

$$Q = m \times c_p \times \Delta T$$

$$Q = 25 \times 4.18 \times 7$$

$$Q = 731.5 \text{ J}$$

Note the mass is the mass of the copper sulfate solution only. Do not include mass of zinc powder.

Step 2: calculate the number of moles of the reactant not in excess.

$$\begin{aligned} \text{moles of CuSO}_4 &= \text{conc} \times \text{vol} \\ &= 0.2 \times 25/1000 \\ &= 0.005 \text{ mol} \end{aligned}$$

If you are not told what is in excess, then you need to work out the moles of both reactants and work out using the balanced equation which one is in excess.

Step 3: calculate the enthalpy change per mole,  $\Delta H$  (the enthalpy change of reaction)

$$\begin{aligned} \Delta H &= Q / \text{no of moles} \\ &= 731.5 / 0.005 \\ &= 146300 \text{ J mol}^{-1} \\ &= 146 \text{ kJ mol}^{-1} \text{ to 3 sf} \end{aligned}$$

Finally add in the sign to represent the energy change: if temp increases the reaction is exothermic and is given a minus sign e.g.  $-146 \text{ kJ mol}^{-1}$

Remember in these questions: **sign, unit**

**Example 4.**  $25.0 \text{ cm}^3$  of  $2.0 \text{ mol dm}^{-3}$  HCl was neutralised by  $25.0 \text{ cm}^3$  of  $2.0 \text{ mol dm}^{-3}$  NaOH. The temperature increased  $13.5^\circ\text{C}$ . Calculate the enthalpy change per mole of HCl?

Step 1: Calculate the energy change for the amount of reactants in the test tube.

$$Q = m \times c_p \times \Delta T$$

$$Q = 50 \times 4.18 \times 13.5$$

$$Q = 2821.5 \text{ J}$$

Note the mass equals the mass of acid + the mass of alkali, as they are both solutions.

Step 2: calculate the number of moles of the HCl.

$$\begin{aligned} \text{moles of HCl} &= \text{conc} \times \text{vol} \\ &= 2 \times 25/1000 \\ &= 0.05 \text{ mol} \end{aligned}$$


Step 3: calculate  $\Delta H$  the enthalpy change per mole which can be called the enthalpy change of neutralisation

$$\begin{aligned} \Delta H &= Q / \text{no of moles} \\ &= 2821.5 / 0.05 \\ &= 56430 \text{ J mol}^{-1} \\ &= -56.4 \text{ kJ mol}^{-1} \text{ to 3 sf} \end{aligned}$$

Exothermic and so is given a minus sign

Remember in these questions: **sign, unit,**

Often Hess's law cycles are used to measure the enthalpy change for a reaction that cannot be measured directly by experiments. Instead alternative reactions are carried out that can be measured experimentally.



This Hess's law is used to work out the enthalpy change to form a hydrated salt from an anhydrous salt.

This cannot be done experimentally because it is impossible to add the exact amount of water without the solid dissolving and it is not easy to measure the temperature change of a solid.

Instead both salts are dissolved in excess water to form a solution of copper sulfate. The temperature changes can be measured for these reactions.

#### Detailed method for measuring enthalpy change of solution of anhydrous copper(II) sulfate

1. Weigh out between 3.90 g and 4.10 g of anhydrous copper(II) sulfate in a dry weighing bottle. The precise mass should be recorded.
2. Using a volumetric pipette, place 25 cm<sup>3</sup> of deionised water into a polystyrene cup and record its temperature at the beginning (t=0), start the timer and then record the temperature again every minute, stirring the liquid continuously.
3. At the fourth minute, add the powdered anhydrous copper(II) sulfate rapidly to the water in the polystyrene cup and continue to stir, but do not record the temperature.
4. Reweigh the empty weighing bottle
5. At the fifth minute and for every minute up to 15 minutes, stir and record the temperature of the solution in the polystyrene cup.
6. Plot a graph of temperature (on the y-axis) against time. Draw two separate best fit lines; one, which joins the points before the addition, and one, which joins the points after the addition, extrapolating both lines to the fourth minute.
7. Use your graph to determine the temperature change at the fourth minute, which theoretically should have occurred immediately on addition of the solid.
8. Using  $q = m \times c_p \times \Delta T$  calculate energy change  
 $= 20 \times 4.18 \times \Delta T$
9. Calculate  $\Delta H_{\text{solution}}$  by dividing q by number of moles of anhydrous copper(II) sulfate in mass added

The above method is then repeated using hydrated copper sulfate. The two  $\Delta H_{\text{solution}}$  can then be used to calculate the  $\Delta H$  for the enthalpy change of forming a hydrated salt as in the example above

## Measuring Enthalpies of Combustion using Flame Calorimetry


Enthalpies of combustion can be calculated by using calorimetry. Generally the fuel is burnt and the flame is used to heat up water in a metal cup.

Need to measure

- mass of spirit burner before and after
- Temperature change of water
- Volume of water in cup

### Errors in this method

- energy losses from calorimeter
- Incomplete combustion of fuel
- Incomplete transfer of energy
- Evaporation of fuel after weighing
- Heat capacity of calorimeter not included
- Measurements not carried out under standard conditions as H<sub>2</sub>O is gas, not liquid, in this experiment



**Example 5.** Calculate the enthalpy change of combustion for the reaction where 0.65g of propan-1-ol was completely combusted and used to heat up 150g of water from 20.1 to 45.5°C

Step 1: Calculate the energy change used to heat up the water.

$$Q = m \times c_p \times \Delta T$$

$$Q = 150 \times 4.18 \times 25.4$$

$$Q = 15925.8 \text{ J}$$

Note the mass in this equation is the mass of water in the calorimeter and not the alcohol

Step 2 : calculate the number of moles of alcohol combusted.

$$\begin{aligned} \text{moles of propan-1-ol} &= \text{mass} / M_r \\ &= 0.65 / 60 \\ &= 0.01083 \text{ mol} \end{aligned}$$

Step 3 : calculate the enthalpy change per mole which is called  $\Delta_c H$  (the enthalpy change of combustion)

$$\begin{aligned} \Delta H &= Q / \text{no of moles} \\ &= 15925.8 / 0.01083 \\ &= 1470073 \text{ J mol}^{-1} \\ &= 1470 \text{ kJ mol}^{-1} \text{ to 3 sf} \end{aligned}$$

Finally add in the sign to represent the energy change: if temp increases the reaction is exothermic and is given a minus sign eg **-1470 kJ mol<sup>-1</sup>**


Remember in these questions: **sign, unit**

## Core practical 3+4. Make up a volumetric solution and carry out a simple acid–base titration

### Making a solution

- Weigh the sample bottle containing the required mass of solid on a 2 dp balance
  - Transfer to beaker
  - Reweigh empty sample bottle
  - Record the difference in mass
- 
- Add 100cm<sup>3</sup> of distilled water to the beaker. Use a glass rod to stir to help dissolve the solid.
- 
- Sometimes the substance may not dissolve well in cold water so the beaker and its contents could be heated gently until all the solid had dissolved.
- 
- Pour solution into a 250cm<sup>3</sup> graduated flask via a funnel.
  - Rinse beaker and funnel and add washings from the beaker and glass rod to the volumetric flask.
  - make up to the mark with distilled water using a dropping pipette for last few drops.
  - Invert flask several times to ensure uniform solution.

Alternatively the known mass of solid in the weighing bottle could be transferred to beaker, washed and washings added to the beaker.



Remember to fill so the bottom of the meniscus sits on the line on the neck of the flask. With dark liquids like potassium manganate it can be difficult to see the meniscus.

Shake the volumetric flask thoroughly to ensure a uniform concentration



Graduated/volumetric flask

A graduated flask has one mark on the neck which the level to fill to get the accurate volume. Do not heat or put hot solutions in the volumetric flask because the heat would cause the flask to expand and the volume would then be incorrect.

#### Measuring mass accurately:

In many experiments the best method for measuring mass is

1. Measure mass on 2 or 3d.p. balance of a weighing bottle with the required quantity of solid in it
2. Empty mass into reaction vessel/flask
3. Reweigh the now empty weighing bottle
4. Subtract the mass of the empty weighing bottle from the first reading to give exact of mass actually added.

### Dilutions

#### Diluting a solution

- Pipette 25cm<sup>3</sup> of original solution into a 250cm<sup>3</sup> volumetric flask
- make up to the mark with distilled water using a dropping pipette for last few drops.
- Invert flask several times to ensure uniform solution.

Using a volumetric pipette is more accurate than a measuring cylinder because it has a smaller uncertainty

Use a teat pipette to make up to the mark in volumetric flask to ensure volume of solution accurately measured and one doesn't go over the line



## Titration Core Practical 3+4: Make up a volumetric solution and carry out a simple acid–base titration

### General Method

- rinse equipment (burette with acid, pipette with alkali, conical flask with distilled water)
- pipette 25 cm<sup>3</sup> of alkali into conical flask
- touch surface of alkali with pipette ( to ensure correct amount is added)
- adds acid solution from burette
- make sure the jet space in the burette is filled with acid
- add a few drops of indicator and refer to colour change at end point
- phenolphthalein [pink (alkali) to colourless (acid): end point pink colour just disappears] [use if NaOH is used]
- methyl orange [yellow (alkali) to red (acid): end point orange] [use if HCl is used]
- use a white tile underneath the flask to help observe the colour change
- add acid to alkali whilst **swirling the mixture** and **add acid drop wise at end point**
- note burette reading before and after addition of acid
- repeats titration until **at least 2 concordant results** are obtained- two readings within 0.1 of each other

Titration is done often to find out the concentration of one substance by reacting it with another substance of known concentration.


They are often done with neutralisation reactions, but can be done with redox reactions.

One substance (generally the one we don't know the concentration) is put in the conical flask. It is measured using a volumetric pipette.

The other substance is placed in the burette

However, the standard phrase: **titrate solution A with solution B** means that A should be in the conical flask and B should be in the burette.

A conical flask is used in preference to a beaker because it is easier to swirl the mixture in a conical flask without spilling the contents.




### Detailed Method for Titration

#### Using the pipette

- rinse pipette with substance to go in it (often alkali).
- pipette 25 cm<sup>3</sup> of solution A into conical flask. The volumetric pipette will have a mark on its neck to show the level to fill to. The bottom of the meniscus should sit on this line.
- touch surface of solution with pipette ( to ensure correct amount is added). A small amount of solution will be left in the pipette at this stage. The calibration of the pipette will take into account this effect. It should not be forced out.

Make sure bottom of meniscus is on line on neck of pipette



## Using the burette

The burette should be rinsed out with substance that will be put in it. If it is not rinsed out the acid or alkali added may be diluted by residual water in the burette or may react with substances left from a previous titration. This would lead to the concentration of the substance being lowered and a larger titre being delivered.

Don't leave the funnel in the burette because small drops of liquid may fall from the funnel during the titration leading to a false burette reading (would give a lower titre volume)

**make sure the jet space** in the burette is **filled** with the solution and air bubbles are removed.


If the jet space in the burette is not filled properly prior to commencing the titration it will lead to errors if it then fills during the titration, leading to a larger than expected titre reading.



Read the bottom of the meniscus on the burette

This is reading  $9.00\text{cm}^3$

Even though a burette has marking reading to  $0.1\text{cm}^3$ , the burette readings should always be given to 2dp either ending in  $0.00$  or  $0.05$ .  $0.05\text{cm}^3$  is the volume of 1 drop of solution delivered from a burette and so this is the smallest difference in readings that can be measured. If the bottom of the meniscus sits on a line it should end with a  $0.00$  as in the above example  $9.00\text{cm}^3$ . If the meniscus sits between two lines it should end  $0.05$ . e.g. if the bottom of the meniscus sits between the lines marked 9.1 and 9.2, you should record  $9.15$



## Adding indicator

Add **a few drops of indicator** and refer to colour change at end point

Indicators are generally weak acids so only add a few drops of them. If too much is added they will affect the titration result

### phenolphthalein

If acid is added from the burette the colour change would be pink (alkali) to colourless (acid): end point pink colour just disappears [use with titrations using strong alkalis e.g. NaOH ]



phenolphthalein  
Alkali colour



phenolphthalein acid  
colour

### Methyl orange

Methyl orange is a suitable indicator for neutralisation reactions where strong acids are used. It is red in acid and yellow in alkali. It is orange at the end point.



Methyl orange  
Alkali colour



Methyl orange  
end point



Methyl orange  
acid colour

Use a white tile underneath the flask to help observe the colour change

Add solution from burette whilst **swirling the mixture** and **add drop-wise at end point**

Distilled water can be added to the conical flask during a titration to wash the sides of the flask so that all the acid on the side is washed into the reaction mixture to react with the alkali.  
It does not affect the titration reading as water does not react with the reagents or change the number of moles of acid added.

**note burette reading** before and after addition of solution  
**repeats titration until at least 2 concordant results** are obtained- two readings within 0.1 of each other

A single titration could be flawed. Repeating allows for anomalous titres to be spotted and discounted

#### Recording results

- Results should be clearly recorded in a table
- Result should be recorded in full (i.e. both initial and final readings)
- Record titre volumes to 2dp (0.05 cm<sup>3</sup>)

Titration number	1	2	3
Initial burette reading (cm <sup>3</sup> )	0.50	2.50	1.55
Final burette reading (cm <sup>3</sup> )	24.50	27.00	25.95
Titre (cm <sup>3</sup> )	24.00	24.50	24.40

#### Safety precautions

**Acids and alkalis are corrosive  
(at low concentrations acids are irritants)**

**Wear eye protection and gloves**

If spilled immediately wash affected parts after spillage

If substance is unknown treat it as potentially toxic and wear gloves.

#### Testing batches

In quality control it will be necessary to do titrations/testing on several samples as the amount/concentration of the chemical being tested may vary between samples.

Only distilled water should be used to wash out conical flasks between titrations because it does not add any extra moles of reagents






If **2 or 3 values are within 0.10cm<sup>3</sup>** and therefore **concordant** or close then we can say results are accurate and **repeatable** and **the titration technique is good and consistent**

#### Working out average titre results

Only make an average of the concordant titre results

$$\text{Average titre} = (24.50 + 24.40) / 2 = 24.45$$

#### Common Titration Equations



#### Titration mixtures

If titrating a mixture to work out the concentration of an active ingredient it is necessary to consider if the mixture contains other substances that have acid base properties.

If they don't have acid base properties we can titrate with confidence.

## Uncertainty

### Readings and Measurements

#### Readings

the values found from a single judgement when using a piece of equipment

#### Measurements

the values taken as the difference between the judgements of two values (e.g. using a burette in a titration)

In general, if uncertainty is not indicated on apparatus, the following assumptions are made:

For an analogue scale-

The uncertainty of a reading (one judgement) is at least  $\pm 0.5$  of the smallest scale reading.

The uncertainty of a measurement (two judgements) is at least  $\pm 1$  of the smallest scale reading.

- If the apparatus has a digital scale, the uncertainty is  $\pm$  the resolution of the apparatus in each measurement

#### Calculating Apparatus Uncertainties

Each type of apparatus has a sensitivity uncertainty

- balance  $\pm 0.001$  g (if using a 3 d.p. balance)
- volumetric flask  $\pm 0.1$  cm<sup>3</sup>
- 25 cm<sup>3</sup> pipette  $\pm 0.1$  cm<sup>3</sup>
- burette (start & end readings and end point)  $\pm 0.15$  cm<sup>3</sup>

Calculate the percentage error for each piece of equipment used by

$$\% \text{ uncertainty} = \pm \frac{\text{uncertainty}}{\text{Measurement made on apparatus}} \times 100$$

e.g. for burette

$$\% \text{ uncertainty} = 0.15 / \text{average titre result} \times 100$$

To calculate the maximum **total** percentage apparatus uncertainty in the final result add all the individual equipment uncertainties together.

#### Uncertainty of a measurement using a burette.

If the burette used in the titration had an uncertainty for each reading of  $\pm 0.05$  cm<sup>3</sup> then during a titration two readings would be taken so the uncertainty on the titre volume would be  $\pm 0.10$  cm<sup>3</sup>.

#### Reducing uncertainties in a titration

Replacing measuring cylinders with pipettes or burettes which have lower apparatus uncertainty will lower the % uncertainty.

To reduce the % uncertainty in a burette reading it is necessary to make the titre a larger volume. This could be done by: increasing the volume and concentration of the substance in the conical flask or by decreasing the concentration of the substance in the burette.

To decrease the apparatus uncertainties you can either decrease the sensitivity uncertainty by using apparatus with a greater resolution (finer scale divisions) or you can increase the size of the measurement made.

#### Reducing uncertainties in measuring mass

Using a balance that measures to more decimal places or using a larger mass will reduce the % uncertainty in weighing a solid.

Weighing sample before and after addition and then calculating difference will ensure a more accurate measurement of the mass added.

If looking at a series of measurements in an investigation, the experiments with the smallest readings will have the highest experimental uncertainties.

#### Calculating the percentage difference between the actual value and the calculated value

If we calculated an Mr of 203 and the real value is 214, then the calculation is as follows:

Calculate difference  $214 - 203 = 11$

$$\% = 11 / 214 \times 100$$

$$= 5.41\%$$

If the %**uncertainty** due to the apparatus  $<$  percentage difference between the actual value and the calculated value then there is a discrepancy in the result due to other errors.

If the %**uncertainty** due to the apparatus  $>$  percentage difference between the actual value and the calculated value then there is no discrepancy and all errors in the results can be explained by the sensitivity of the equipment.

## Core practical 5: Investigation of the rates of hydrolysis of some halogenoalkanes

### Testing for halogenoalkanes method

- Arrange three test tubes in a row and add three drops of halogenoalkane in the sequence 1-chlorobutane, 1-bromobutane, 1-iodobutane.
- Add 4 cm<sup>3</sup> of 0.02 M silver nitrate to each halogenoalkane.
- Without delay, put all three test tubes simultaneously in a hot water bath.
- Note the order in which precipitates appear

### Comparing the rate of hydrolysis of halogenoalkanes reaction

**Hydrolysis** is defined as the splitting of a molecule ( in this case a haloalkane) by a reaction with water





**Aqueous silver nitrate** is added to a haloalkane and the halide leaving group combines with a silver ion to form a **SILVER HALIDE PRECIPITATE**.

The precipitate only forms when the halide ion has left the haloalkane and so the rate of formation of the precipitate can be used to compare the reactivity of the different haloalkanes.

The **quicker** the **precipitate** is formed, the **faster** the **substitution** reaction and the **more reactive** the haloalkane

The rate of these substitution reactions depends on the strength of the C-X bond . The weaker the bond, the easier it is to break and the faster the reaction.

**Water** is a **poor nucleophile** but it can react **slowly** with haloalkanes in a substitution reaction



The iodoalkane forms a precipitate with the silver nitrate first as the C-I bond is weakest and so it hydrolyses the quickest

AgI<sub>(s)</sub> - yellow precipitate

AgBr<sub>(s)</sub> - cream precipitate

AgCl<sub>(s)</sub> - white precipitate


↑ forms faster

## Purifying an organic liquid

### General method

- Put the distillate of impure product into a separating funnel
- wash product by adding either
  - sodium hydrogencarbonate solution, shaking and releasing the pressure from CO<sub>2</sub> produced.
  - Saturated sodium chloride solution
- Allow the layers to separate in the funnel, and then run and discard the aqueous layer.
- Run the organic layer into a clean, dry conical flask and add three spatula loads of drying agent (e.g. anhydrous sodium sulphate, calcium chloride) to dry the organic liquid. When dry the organic liquid should appear clear.
- Carefully decant the liquid into the distillation flask
- Distill to collect pure product

Separating funnel



Sodium hydrogencarbonate will neutralise any remaining reactant acid.

Sodium chloride will help separate the organic layer from the aqueous layer

The layer with lower density will be the upper layer. This is usually the organic layer

The drying agent should

- be insoluble in the organic liquid
- not react with the organic liquid

Decant means carefully pour off organic liquid leaving the drying agent in the conical flask


### Distillation

In general used as separation technique to separate an organic product from its reacting mixture. Need to collect the distillate of the approximate boiling point range of the desired liquid.

Note the bulb of the thermometer should be at the T junction connecting to the condenser to measure the correct boiling point

Note the water goes in the bottom of the condenser to go against gravity. This allows more efficient cooling and prevents back flow of water.

Electric heaters are often used to heat organic chemicals. This is because organic chemicals are normally highly flammable and could set on fire with a naked flame.




## CORE PRACTICAL 6: Chlorination of 2-methylpropan-2-ol using concentrated hydrochloric acid

### Detailed method for preparing and purifying a halogenoalkane from an alcohol

1. Measure 8 cm<sup>3</sup> of 2-methylpropan-2-ol in a measuring cylinder and measure its mass.
2. Pour the 2-methylpropan-2-ol into a separating funnel, and reweigh the measuring cylinder to find the mass of the 2-methylpropan-2-ol used.
3. In a fume cupboard, add 20 cm<sup>3</sup> of concentrated hydrochloric acid to the separating funnel, in portions of 3cm<sup>3</sup>. After each portion, stopper the flask and invert it several times. Open the tap when doing this to release the pressure.
4. Allow the separating funnel to stand in the fume cupboard for about 20 minutes. Gently shake it at intervals.
5. After 20 minutes, allow the layers to separate in the funnel. Open the tap and remove the lower aqueous layer. Dispose of this layer.
6. Add sodium hydrogencarbonate solution in 2 cm<sup>3</sup> portions to the separating funnel. This neutralises any remaining acid. Shake the funnel after each addition, and release the pressure. Continue until no more bubbles of CO<sub>2</sub> are seen.
7. Allow the layers to separate in the funnel. Again remove and pour away the lower aqueous layer. Run off the organic layer into a clean conical flask and add two spatulas of anhydrous sodium sulfate. Stopper the flask, shake the contents and allow this to stand until the liquid becomes clear. This step dries the organic liquid.
8. Decant the liquid into a weighed clean distillation flask.
9. Distil the liquid by holding a 250ml beaker half-full of boiled water around the flask using standard distillation set up. Collect the liquid that distils in the range 47-53°C.
10. Measure the mass of the 2-chloro-2-methylpropane collected.

### Detailed method for preparing and purifying Cyclohexene from cyclohexanol




- a) Pour about 20 cm<sup>3</sup> of cyclohexanol into a weighed 50 cm<sup>3</sup> pear-shaped flask. Reweigh the flask and record the mass of cyclohexanol.
- b) Using a plastic graduated dropping pipette, carefully and with frequent shaking, add to the flask approximately 8 cm<sup>3</sup> of concentrated phosphoric acid.
- c) Add a few anti-bumping granules to the flask and assemble the distillation apparatus, so that the contents of the flask may be distilled. Heat the flask gently, distilling over any liquid which boils below 100 °C.
- d) Pour the distillate into a separating funnel and add 50 cm<sup>3</sup> of saturated sodium chloride solution. Shake the mixture and allow the two layers to separate.
- e) run off the lower layer into a beaker and then transfer the upper layer, which contains the crude cyclohexene, into a small conical flask.
- f) Add a few lumps of anhydrous calcium chloride or anhydrous sodium sulfate(VI) or anhydrous magnesium sulfate to the crude cyclohexene to remove water. Stopper the flask, shake the contents and allow this to stand until the liquid becomes clear.
- g) Decant the liquid into a clean, dry, weighed sample container.
- h) Reweigh the container, calculate the mass of dry cyclohexene produced

### Detailed Method for Preparing and Purifying an Ester

Propyl ethanoate can be made in the laboratory from propan-1-ol and ethanoic acid.

The equation for the reaction is



#### Procedure

1. Propan-1-ol (50 cm<sup>3</sup>) and ethanoic acid (50 cm<sup>3</sup>) are mixed thoroughly in a 250 cm<sup>3</sup> round-bottomed flask.
2. **Concentrated sulfuric acid** (10 cm<sup>3</sup>) is added **drop by drop to the mixture**, keeping the contents of the flask well-shaken and **cooled** in an ice-water bath.
3. When the acid has all been added, a **reflux condenser** is fitted to the flask and the mixture gently boiled over an **electric heating mantle** for about **30 minutes**.
4. The mixture is cooled, and the apparatus rearranged for distillation. The crude ester (about 60 cm<sup>3</sup>) is distilled off.
5. The distillate is placed in a separating funnel and shaken with about half its volume of **30% sodium carbonate** solution, with the **pressure being released at intervals**. The **lower aqueous layer is then discarded**.
6. The crude ester is shaken in a separating funnel with about half its volume of 50% calcium chloride solution, which removes unreacted alcohol. The lower layer is discarded.
7. The ester is run into a clean, dry flask containing **some anhydrous calcium chloride** and swirled.
8. The ester is filtered into a clean, dry flask, with a **few anti-bumping granules**, and distilled. The fraction boiling between 100°C and 103°C is collected.

**Sulfuric acid** is a catalyst

Adding conc H<sub>2</sub>SO<sub>4</sub> is an exothermic reaction- to prevent uncontrolled boiling over add **drop by drop** and **cool**

In **reflux** the reactant vapours of volatile compound are condensed and returned to the reaction mixture. The reaction is slow so it is heated for **30 minutes** **The electric heating mantle** allows for controlled heating and stops flammable vapour lighting

**Sodium carbonate** reacts with unreacted acid and remaining catalyst still present after distillation.

The reaction produces CO<sub>2</sub> so the **pressure** of gas needs to be **released**.

The upper layer is organic because it has a lower density than water

**Calcium chloride** is a drying agent. The liquid will appear clear when dry.

**Anti-bumping granules** are added to prevent vigorous, uneven boiling by making small bubbles form instead of large bubbles

### Measuring boiling point

Purity of liquid can be determined by measuring a boiling point. This can be done in a distillation set up or by simply boiling a tube of the sample in an heating oil bath. If the liquid is pure it will have the boiling point referred to in data books. If impure the boiling point tends to be higher than the pure liquid

Pressure should be noted as changing pressure can change the boiling point of a liquid

Measuring boiling point is not the most accurate method of identifying a substance as several substances may have the same boiling point.

To get a correct measure of boiling point the thermometer should be above the level of the surface of the boiling liquid and be measuring the temperature of the saturated vapour.




## Core practical 7. Oxidation of an alcohol

### Partial Oxidation of Primary Alcohols

**Reaction:** primary alcohol  $\rightarrow$  aldehyde

**Reagent:** potassium dichromate (VI) solution and dilute sulfuric acid.

**Conditions:** (use a limited amount of dichromate) warm gently and **distil** out the aldehyde as it forms:




Ethanal

**Observation:** the orange dichromate ion ( $\text{Cr}_2\text{O}_7^{2-}$ ) reduces to the green  $\text{Cr}^{3+}$  ion

#### SAFETY

You must wear gloves when handling solid potassium dichromate(VI) since it is highly toxic and a category 2 carcinogen; it is also an irritant. Avoid inhaling any dust.

Concentrated sulfuric acid is corrosive.



#### Detailed Method: The partial oxidation of propan-1-ol

This experiment uses a limited quantity of oxidising agent (0.01 mol) and the product is distilled from the reaction mixture immediately it is formed. In this way we hope to achieve a partial oxidation of propan-1-ol.

- Place about 10 cm<sup>3</sup> of dilute sulfuric acid in a flask and add about 3g of potassium dichromate(VI) and 2 or 3 anti-bumping granules. Shake the contents of the flask until solution is complete (do not warm).
- Add 1.5 cm<sup>3</sup> of propan-1-ol in drops from a dropping pipette, shaking the flask so as to mix the contents, and then assemble distillation apparatus as shown below
- Gently heat and slowly distil 2 cm<sup>3</sup> of liquid into a test tube, taking care that none of the reaction mixture splashes over.


#### Distillation

In general used as separation technique to separate an organic product from its reacting mixture. Need to collect the distillate of the approximate boiling point range of the desired liquid.

Note the bulb of the thermometer should be at the T junction connecting to the condenser to measure the correct boiling point

Note the water goes in the bottom of the condenser to go against gravity. This allows more efficient cooling and prevents back flow of water.

Electric heaters are often used to heat organic chemicals. This is because organic chemicals are normally highly flammable and could set on fire with a naked flame.




It's important to be able to draw and label this apparatus accurately. Don't draw lines between flask, adaptor and condenser.



### Fractional Distillation: In the laboratory

- Heat the flask, with a Bunsen burner or electric mantle
- This causes vapours of all the components in the mixture to be produced.
- Vapours pass up the fractionating column.
- The vapour of the substance with the lower boiling point reaches the top of the fractionating column first.
- The thermometer should be at or below the boiling point of the most volatile substance.
- The vapours with higher boiling points condense back into the flask.
- Only the most volatile vapour passes into the condenser.
- The condenser cools the vapours and condenses to a liquid and is collected.

Fractional distillation is used to separate liquids with different boiling points



## CORE PRACTICAL 8 + 15: Analysis of some inorganic unknowns

### Testing for cations

#### Flame tests

##### Method

Use a **nichrome wire** ( nichrome is an unreactive metal and will not give out any flame colour)  
Clean the wire by dipping in **concentrated hydrochloric acid** and then heating in Bunsen flame  
If the sample is not powdered then grind it up.  
Dip wire in solid and put in Bunsen flame and observe flame

Lithium : Scarlet red  
Sodium : Yellow  
Potassium : lilac  
Rubidium : red  
Caesium: blue  
Magnesium: no flame colour (energy emitted of a wavelength outside visible spectrum)  
Calcium: brick red  
Strontium: red  
Barium: apple green

##### Method: adding dilute sodium hydroxide

- Place about 10 drops of 0.1 mol dm<sup>-3</sup> metal ion solution in a test tube.
- Add about 10 drops of 0.6 mol dm<sup>-3</sup> sodium hydroxide solution, mixing well.
- Continue to add sodium hydroxide solution, dropwise with gentle shaking, until in excess


This test can be used on group 2 metal ions and transition metal ions.

#### Results for Group 2

The results in this test are an application of the trend that group II hydroxides become **more soluble** down the group.

**Magnesium hydroxide** is classed as insoluble in water and will appear as a **white precipitate**.


**Simplest Ionic Equation** for formation of Mg(OH)<sub>2</sub> (s)



A suspension of magnesium hydroxide in water will appear slightly alkaline (pH 9) so some hydroxide ions must therefore have been produced by a very slight dissolving.

**Calcium hydroxide** is classed as partially soluble in water and will appear as a **white precipitate** (it may need more sodium hydroxide to be added before it appears compared to a magnesium solution.)

**Simplest Ionic Equation** for formation of Ca(OH)<sub>2</sub> (s)









A suspension of calcium hydroxide in water will appear more alkaline (pH 11) than magnesium hydroxide as it is more soluble so there will be more hydroxide ions present in solution.

**Strontium** and **barium** salts will not form a hydroxide precipitate on addition of sodium hydroxide due to their high solubility. The solutions will be highly alkaline

#### Results for transition metals

Copper(II) solutions form a blue ppt,  
Cobalt(II) solutions form a blue ppt,  
iron (II) solutions form a green ppt  
iron (III) solutions form a brown ppt  
Chromium (III) solutions form a green ppt which dissolves in excess to form a green solution




#### Testing for Ammonium ions (NH<sub>4</sub><sup>+</sup>)


- Place about 10 drops of 0.1 mol dm<sup>-3</sup> ammonium chloride in a test tube.
- Add about 10 drops of 0.4 mol dm<sup>-3</sup> sodium hydroxide solution. Shake the mixture.
- Warm the mixture in the test tube gently using a water bath.
- Test the fumes released from the mixture by holding a piece of damp red litmus paper in the mouth of the test tube.

Results: alkaline ammonia gas is released which turns the red litmus paper blue

## 2+ Ion Summary



## 3+ Ion Summary



## CORE PRACTICAL 8+15: Analysis of some inorganic unknowns


### Testing for anions: – Group 7 (halide ions), OH<sup>-</sup>, CO<sub>3</sub><sup>2-</sup>, SO<sub>4</sub><sup>2-</sup>

#### Testing for presence of a sulfate ion

**BaCl<sub>2</sub> solution acidified with hydrochloric acid** is used as a reagent to test for sulphate ions.

If acidified **Barium Chloride** is added to a solution that contains sulfate ions a **white precipitate** of Barium Sulfate forms.


Simplest ionic equation



Other anions should give a negative result which is no precipitate forming.



The hydrochloric acid is needed to react with carbonate impurities that are often found in salts which would form a white Barium carbonate precipitate and so give a false result. You could not use sulphuric acid because it contains sulphate ions and so would give a false positive result.




Fizzing due to CO<sub>2</sub> would be observed if a carbonate was present.

#### Testing for presence of halide ions with silver nitrate.


This reaction is used as a test to identify which halide ion is present. The test solution is made acidic with **nitric acid**, and then **Silver nitrate solution** is added dropwise.

Fluorides produce no precipitate


Chlorides produce a **white precipitate**




Bromides produce a **cream precipitate**



Iodides produce a **pale yellow precipitate**




The role of nitric acid is to react with any carbonates present to prevent formation of the precipitate Ag<sub>2</sub>CO<sub>3</sub>. This would mask the desired observations




The silver halide precipitates can be treated with ammonia solution to help differentiate between them if the colours look similar:

**Silver chloride** dissolves in **dilute ammonia** to form a complex ion



Colourless solution

**Silver bromide** dissolves in **concentrated ammonia** to form a complex ion




Colourless solution

**Silver iodide** does not react with ammonia – it is too insoluble.

#### Testing for presence of carbonate ions

Add any dilute acid and observe effervescence.

Bubble gas through limewater to test for CO<sub>2</sub> – will turn limewater cloudy



Fizzing due to CO<sub>2</sub> would be observed if a carbonate was present

#### Testing for presence of a hydroxide ions

Alkaline hydroxide ions will turn red litmus paper blue

## Reactions of halide salts with concentrated sulfuric acid.

The halides show increasing power as reducing agents as one goes down the group. This can be clearly demonstrated in the various reactions of the solid halides with concentrated sulfuric acid.


**Know the equations and observations of these reactions very well.**

### Explanation of differing reducing power of halides


A reducing agent donates electrons.  
The reducing power of the halides increases down group 7  
They have a greater tendency to donate electrons.  
This is because as the ions get bigger it is easier for the outer electrons to be given away as the pull from the nucleus on them becomes smaller.

### Fluoride and Chloride

F<sup>-</sup> and Cl<sup>-</sup> ions are not strong enough reducing agents to reduce the S in H<sub>2</sub>SO<sub>4</sub>. No redox reactions occur. Only acid-base reactions occur.



Observations: White steamy fumes of HF are evolved.







Observations: White steamy fumes of HCl are evolved.

These are acid-base reactions and not redox reactions. H<sub>2</sub>SO<sub>4</sub> plays the role of an acid (proton donor).

### Bromide

Br<sup>-</sup> ions are stronger reducing agents than Cl<sup>-</sup> and F<sup>-</sup> and after the initial acid-base reaction reduce the sulfur in H<sub>2</sub>SO<sub>4</sub> from +6 to +4 in SO<sub>2</sub>



**Observations:** White steamy fumes of HBr are evolved.









Red fumes of Bromine are also evolved and a colourless, acidic gas SO<sub>2</sub>

Reduction product = sulfur dioxide

Note the H<sub>2</sub>SO<sub>4</sub> plays the role of acid in the first step producing HBr and then acts as an oxidising agent in the second redox step.

### Iodide

I<sup>-</sup> ions are the strongest halide reducing agents. They can reduce the sulfur from +6 in H<sub>2</sub>SO<sub>4</sub> to +4 in SO<sub>2</sub>, to 0 in S and -2 in H<sub>2</sub>S.



### Observations:

White **steamy fumes** of HI are evolved.

**Black solid** and **purple fumes** of Iodine are also evolved

A **colourless**, acidic gas SO<sub>2</sub>

A **yellow solid** of sulfur

**H<sub>2</sub>S** (Hydrogen sulfide), a **gas** with a **bad egg smell**.

Reduction products = sulfur dioxide, sulfur and hydrogen sulfide

Note the H<sub>2</sub>SO<sub>4</sub> plays the role of acid in the first step producing HI and then acts as an oxidising agent in the three redox steps

*Often in exam questions these redox reactions are worked out after first making the half-equations*

## More on Insoluble salts and Precipitation reactions


Insoluble salts can be made by mixing appropriate solutions of ions so that a **precipitate** is formed  
 Barium nitrate (aq) + sodium sulfate (aq) → **barium sulfate (s)** + sodium nitrate (aq)  
 These are called **precipitation** reactions. A **precipitate is a solid**

There are some common rules for solubility of salts. No syllabus requires these to be learnt but a good chemist does know them.


Soluble salts	Insoluble salts
All sodium, potassium and ammonium salts	
All nitrates	
Most chlorides, bromides, iodides	Silver, lead chlorides, bromides iodides
Most sulfates	Lead, strontium and barium sulfate
Sodium, potassium and ammonium carbonates	Most other carbonates
Sodium, potassium and ammonium hydroxides	Most other hydroxides

When making an insoluble salt, normally the salt would be removed by **filtration**, washed with **distilled water to remove soluble impurities** and then **dried on filter paper**

### Filtration



This is gravitational filtration. Use if small amounts of solid are formed.



This is vacuum filtration. The apparatus is connected to a water pump which will produce a vacuum. Use if larger amounts of solid are formed.

For both types of filtration apparatus draw filter paper on the diagram.


### Writing Ionic equations for precipitation reactions

We usually write ionic equations to show precipitation reactions. Ionic equations only show the ions that are reacting and leave out spectator ions.


Spectator ions are ions that are

- Not changing state
- Not changing oxidation number


Take full equation



Separate (aq) solutions into ions



Cancel out spectator ions leaving the simplest ionic equation





## CORE PRACTICAL 8+15: Analysis of some organic unknowns

### Tests for alcohol, aldehyde, alkene and carboxylic acid

#### Functional group test for an Alkene

To 0.5 cm<sup>3</sup> of bromine water in a test tube add a few drops of the unknown and shake.

**Observation:** alkenes should decolourise bromine water

#### Reaction with 2,4-dinitro phenylhydrazine

2,4-DNP reacts with both aldehydes and ketones. The product is an orange precipitate, It can be used as a test for a carbonyl group in a compound.

Use 2,4-DNP to identify if the compound is a carbonyl. Then to differentiate an aldehyde from a ketone use Tollen's reagent.

The melting point of the crystal formed can be used to help identify which carbonyl was used. Take the melting point of orange crystals product from 2,4-DNP. Compare melting point with known values in database

#### Functional group tests for an aldehyde


##### Tollen's Reagent

**Reagent:** Tollen's reagent formed by mixing aqueous ammonia and silver nitrate. The active substance is the complex ion of  $[\text{Ag}(\text{NH}_3)_2]^+$ .

**Conditions:** heat gently


**Reaction:** **aldehydes only** are oxidised by Tollen's reagent into a carboxylic acid and the silver(I) ions are reduced to silver atoms

**Observation:** with aldehydes, a silver mirror forms coating the inside of the test tube. Ketones result in no change.



##### Tollen's reagent method

Place 1 cm<sup>3</sup> of silver nitrate solution in each of two clean boiling tubes. Then add one drop of sodium hydroxide solution to form a precipitate of silver oxide. Add ammonia solution dropwise until a clear, colourless solution is formed. Add a few drops of the unknown and leave in the water bath for a few minutes.




##### Fehling's solution

**Reagent:** Fehling's solution containing blue  $\text{Cu}^{2+}$  ions.

**Conditions:** heat gently

**Reaction:** **aldehydes only** are oxidised by Fehling's Solution into a carboxylic acid and the copper (II) ions are reduced to copper(I) oxide.

**Observation:** **Aldehydes:** Blue  $\text{Cu}^{2+}$  ions in solution change to a red precipitate of  $\text{Cu}_2\text{O}$ . **Ketones do not react**



##### Fehling's solution method

Place 1 cm<sup>3</sup> of Fehling's A into each of two boiling tubes, and then add Fehling's B until the blue precipitate redissolves. Add a few drops of the unknown and leave in the water bath for a few minutes.


## Reaction of carbonyls with iodine in presence of alkali

**Reagents:** Iodine and sodium hydroxide

**Conditions:** warm very gently

The product  $\text{CHI}_3$  is a yellow crystalline precipitate with an antiseptic smell

Only carbonyls with a methyl group next to the  $\text{C}=\text{O}$  bond will do this reaction. Ethanal is the only aldehyde that reacts. More commonly is methyl ketones.




This reaction is called the Iodoform test

## Functional group test for a Carboxylic acid

To  $0.5 \text{ cm}^3$  of your unknown solution in a test tube add a small amount of sodium carbonate solid and observe.

**Result** carboxylic acids will fizz with sodium carbonate due to  $\text{CO}_2$  produced

The presence of a carboxylic acid can be tested by addition of **sodium carbonate**. It will fizz and produce carbon dioxide



## Summary of Identification of Functional Groups by test-tube reactions

Functional group	Reagent	Result
Alkene	Bromine water	Orange colour decolourises
Alcohols + carboxylic acids	$\text{PCl}_5$	Misty fumes of HCl produced
Alcohols, phenols, carboxylic acids	Sodium metal	Efferevesence due to $\text{H}_2$ gas
Carbonyls	2,4,DNP	Orange/red crystals produced
Aldehyde	Fehling's solution	Blue solution to red precipitate
Aldehyde	Tollen's reagent	Silver mirror formed
Carboxylic acid	Sodium carbonate	Effervescence of $\text{CO}_2$ evolved
1° 2° alcohol and aldehyde	Sodium dichromate and sulfuric acid	Orange to green colour change
chloroalkane	Warm with silver nitrate	Slow formation of white precipitate of AgCl
Acyl chloride	Silver nitrate	Vigorous reaction- steamy fumes of HCl- rapid white precipitate of AgCl




# Investigating rates of reaction

## Techniques to investigate rates of reaction

There are several different methods for measuring reactions rates. Some reactions can be measured in several ways


### measurement of the change in volume of a gas

This works if there is a change in the number of moles of gas in the reaction. Using a gas syringe is a common way of following this.




### Measurement of change of mass

This works if there is a gas produced which is allowed to escape. Works better with heavy gases such as  $\text{CO}_2$




### Titrating samples of reaction mixture with acid, alkali, sodium thiosulfate etc


1. Small samples are removed from the reaction mixture
2. quench (which stops the reaction)- can be done by
  - by dilution with water
  - by cooling
  - by adding a reagent that reacts with one of the reactants
3. then titrate with a suitable reagent.



The NaOH could be titrated with an acid




The  $\text{H}^+$  could be titrated with an alkali



The  $\text{I}_2$  could be titrated with sodium thiosulfate

### Colorimetry.


If one of the reactants or products is coloured then colorimetry can be used to measure the change in colour of the reacting mixtures



The  $\text{I}_2$  produced is a brown solution


### Measuring change in electrical conductivity

Can be used if there is a change in the number of ions in the reaction mixture



### Measurement of optical activity.


If there is a change in the optical activity through the reaction this could be followed in a polarimeter



## General: Measuring the rate of reaction: by an continuous monitoring method



When we follow one experiment over time recording the change in concentration it is the continuous rate method.

The gradient represents the rate of reaction. The reaction is fastest at the start where the gradient is steepest. The rate drops as the reactants start to get used up and their concentration drops. The graph will eventually become horizontal and the gradient becomes zero which represents the reaction having stopped.




### Measurement of the change in volume of a gas

This works if there is a change in the number of moles of gas in the reaction. Using a gas syringe is a common way of following this. It works quite well for measuring continuous rate but a typical gas syringe only measures 100ml of gas so you don't want a reaction to produce more than this volume. Quantities of reactants need to be calculated carefully.



The initial rate is the rate at the start of the reaction, where it is fastest. It is obtained by taking the gradient of a continuous monitoring conc vs time graph at time = zero. A measure of initial rate is preferable as we know the concentrations at the start of the reaction




### Typical Method

- Measure 50 cm<sup>3</sup> of the 1.0 mol dm<sup>-3</sup> hydrochloric acid and add to conical flask.
- Set up the gas syringe in the stand
- Weigh 0.20 g of magnesium.
- Add the magnesium ribbon to the conical flask, place the bung firmly into the top of the flask and start the timer.
- Record the volume of hydrogen gas collected every 15 seconds for 3 minutes.

### Large Excess of reactants

In reactions where there are several reactants, if the concentration of one of the reactant is kept in a large excess then that reactant will appear not to affect rate and will be pseudo-zero order. This is because its concentration stays virtually constant and does not affect rate.

### Continuous rate experiments



### Continuous rate data

This is data from one experiment where the concentration of a substance is followed throughout the experiment.


This data is processed by plotting the data and calculating successive half-lives.

The half-life of a first-order reaction is independent of the concentration and is constant



If half-lives **rapidly** increase then the order is 2<sup>nd</sup> order

## CORE PRACTICAL 9a: Rates of reaction: Following the rate of the iodine propanone reaction by a titrimetric method

Propanone reacts with iodine in acidic solution (the acid is a catalyst) as shown in the equation below.



This reaction can be followed by removing small samples from the reaction mixture with a volumetric pipette. The sample is then quenched by adding excess sodium hydrogencarbonate to neutralize acid catalyst which stops the reaction. Then the sample can be titrated with sodium thiosulphate using a starch catalyst



This reaction is zero order with respect to  $\text{I}_2$  but 1<sup>st</sup> order with respect to the propanone and acid catalyst

The rate equation for the reaction is  
Rate =  $k[\text{CH}_3\text{COCH}_3(\text{aq})][\text{H}^+(\text{aq})]$

If there is a zero order reactant there must be at least two steps in the mechanism because the rate determining step will not involve the zero order reactant

The rate determining step of this reaction must therefore contain one propanone molecule and one  $\text{H}^+$  ion forming an intermediate. The iodine will be involved in a subsequent faster step.

### Detailed Method

1. To a beaker add 25 cm<sup>3</sup> of 1 mol dm<sup>-3</sup> aqueous propanone and 25 cm<sup>3</sup> of 1 mol dm<sup>-3</sup> sulfuric acid.
2. Add 50 cm<sup>3</sup> of 0.02 mol dm<sup>-3</sup> iodine solution. Start the clock swirl the beaker well to mix.
3. Using a 10cm<sup>3</sup> pipette, withdraw a sample of the mixture and transfer it to a conical flask.
4. Add a spatula measure of sodium hydrogencarbonate (This stops the reaction). Record the time at which the sodium hydrogencarbonate is added.
5. Titrate the iodine present in the conical flask with 0.01 mol dm<sup>-3</sup> sodium thiosulfate solution. When the colour turns a pale yellow add the starch indicator. The end point is then when the mixture goes from blue to colourless.
- 6 Every 5 minutes withdraw another 10 cm<sup>3</sup> sample and repeat steps 4 and 5

This method allows the order with respect to iodine to be calculated because the propanone and acid are in large excess so their concentrations do not change during the reaction



## Core practical 9b: Rates of reaction: clock reaction

The initial rate can be calculated from taking the gradient of a continuous monitoring conc vs time graph at time = zero

Initial rate can also be calculated from clock reactions where the time taken to reach a fixed concentration is measured.

### A Common Clock Reaction (no need to learn details)

Hydrogen peroxide reacts with iodide ions to form iodine. The thiosulfate ion then immediately reacts with iodine formed in the second reaction as shown below.



When the  $\text{I}_2$  produced has reacted with all of the limited amount of thiosulfate ions present, excess  $\text{I}_2$  remains in solution. Reaction with the starch then suddenly forms a dark blue-black colour.

A series of experiments is carried out, in which the concentration of iodide ions is varied, while keeping the concentrations of all of the other reagents the same. In each experiment the time taken (t) for the reaction mixture to turn blue is measured.

In clock reactions there are often two successive reactions and an end point is achieved when one limited reactant runs out, resulting in a sudden colour change

By repeating the experiment several times, varying the concentration of a reactant e.g.  $\text{I}^-$ , (keeping the other reactants at constant concentration) you can determine the order of reaction with respect to that reactant

The initial rate of the reaction can be represented as  $(1/t)$

### Detailed method

- Put each of the chemicals in the table in separate burettes.
- In each experiment, measure out required volumes of the potassium iodide, sodium thiosulfate, starch and water into a small conical flask from the burettes
- Measure the hydrogen peroxide into a test tube
- Pour the hydrogen peroxide from the test tube into the conical flask and immediately start the timer. Stir the mixture.
- Time until the first hint of blue/ black colour appears

Experiment	Sulfuric acid ( $\text{H}^+$ ) ml	Starch ml	Water ml	Potassium iodide ( $\text{I}^-$ ) ml	Sodium Thiosulfate $\text{S}_2\text{O}_3^{2-}$ ml	Hydrogen peroxide ml
1	25	1	20	5	5	10
2	25	1	15	10	5	10
3	25	1	10	15	5	10
4	25	1	5	20	5	10
5	25	1	0	25	5	10

Normally to work out the rate equation we do a series of experiments where the initial concentrations of reactants are changed (one at a time) and measure the initial rate each time.

### Working out rate order graphically

In an experiment where the concentration of one of the reagents is changed and the reaction rate measured it is possible to calculate the order graphically

Taking rate equation


$$\text{Rate} = k [\text{I}^-]^n$$

Log both sides of equation

$$\text{Log rate} = \text{log } k + n \text{ log } [\text{Y}]$$

$$Y = c + m x$$



A graph of log rate vs  $\text{log } [\text{I}^-]$  will yield a straight line where the gradient is equal to the order n



In this experiment high concentrations with quick times will have the biggest percentage errors.

## Core practical 10: Finding the activation energy of a reaction

**Experiment** Determine the activation energy for the reaction between bromide ions and bromate(V) ions



The Bromine produced in the first reaction reacts with the phenol. When the phenol is used up, the bromine is no longer removed. The bromine then bleaches the methyl red indicator at the 'end of the reaction'

### Detailed Method

1. Pipette 10 cm<sup>3</sup> of phenol solution and 10 cm<sup>3</sup> of bromide/bromate solution into a boiling tube.
2. Add a few drops of methyl red indicator to the mixture.
3. Into a second boiling tube pipette 5 cm<sup>3</sup> of sulfuric acid.
4. Place the two boiling tubes in a water bath at temperature 20°C.
5. When the substances have reached the water temperature, mix the contents into one of the boiling tubes and swirl. Start the stop clock.
6. Place the boiling tube containing the reaction mixture in the water bath.
7. Stop the clock when the methyl red indicator disappears.
8. Repeat the experiment at 30, 40, 50, 60°C

In this experiment rate is **1/time** where the time is the time taken for the indicator to change colour. This is an approximation for initial rate of reaction as it does not include the change in concentration term. We can use this because we can assume the amount of phenol used in each experiment is the same and constant. The change in concentration is therefore the same for each experiment so only the time taken to reach this concentration is relevant.


### Analysis of results to calculate activation energy

The Arrhenius equation can be rearranged

$$\ln k = \text{constant} - E_a/(RT)$$

k is proportional to the rate of reaction so ln k can be replaced by ln(rate)

From plotting a graph of ln(rate) or ln k against 1/T the activation energy can be calculated from measuring the gradient of the line



### Example 9

Temperature T (K)	1/T	time t (s)	1/t	Ln (1/t)
297.3	0.003364	53	0.018868	-3.9703
310.6	0.00322	24	0.041667	-3.1781
317.2	0.003153	16	0.0625	-2.7726
323.9	0.003087	12	0.083333	-2.4849
335.6	0.00298	6	0.166667	-1.7918

$$\text{gradient} = \frac{y_2 - y_1}{x_2 - x_1}$$


The gradient should always be -ve

In above example gradient = -5680

$$\begin{aligned} E_a &= - \text{gradient} \times R \quad (8.31) \\ &= - (-5680) \times 8.31 \\ &= 47200 \text{ J mol}^{-1} \end{aligned}$$

The unit of EA using this equation will be J mol<sup>-1</sup>. Convert into kJ mol<sup>-1</sup> by dividing 1000

$$E_a = +47.2 \text{ kJ mol}^{-1}$$



use a line of best fit  
use all graph paper  
choose points far apart on the graph to calculate the gradient

## Core practical 11: Finding the Ka value for a weak acid – titration curves

### Constructing a pH curve



1. Transfer 25cm<sup>3</sup> of acid to a conical flask with a volumetric pipette
2. Measure initial pH of the acid with a pH meter
3. Add alkali in small amounts (2cm<sup>3</sup>) noting the volume added
4. Stir mixture to equalise the pH
5. Measure and record the pH to 1 d.p.
6. Repeat steps 3-5 but when approaching endpoint add in smaller volumes of alkali
7. Add until alkali in excess

**Calibrate** meter first by measuring known pH of a buffer solution. This is necessary because pH meters can lose accuracy on storage. Most pH probes are calibrated by putting probe in a set buffer (often pH 4) and pressing a calibration button/setting for that pH. Sometimes this is repeated with a second buffer at a different pH

Can also improve accuracy by **maintaining** constant temperature

### Strong acid – Strong base


e.g. HCl and NaOH



### Weak acid – Strong base

e.g. CH<sub>3</sub>CO<sub>2</sub>H and NaOH

At the start the pH rises quickly and then levels off. The flattened part is called the **buffer region** and is formed because a buffer solution is made



### Half neutralisation volume

For weak acids

$$K_a = \frac{[H^+_{(aq)}][A^-_{(aq)}]}{[HA_{(aq)}]}$$

At ½ the neutralisation volume the [HA] = [A<sup>-</sup>]


$$\text{So } K_a = [H^+] \text{ and } pK_a = pH$$

If we know the K<sub>a</sub> we can then work out the pH at ½ V or vice versa.

If a pH curve is plotted then the pH of a weak acid at half neutralisation (½ V) will equal the pK<sub>a</sub>


### Strong acid – Weak base

e.g. HCl and NH<sub>3</sub>



### Weak acid – Weak base

e.g. CH<sub>3</sub>CO<sub>2</sub>H and NH<sub>3</sub>






## Core practical 12: Investigating some electrochemical cells

### Method

- Clean the zinc and copper foils with emery before use. Degrease the metal using some cotton wool and propanone.
- Place the copper strip into a 100 cm<sup>3</sup> beaker with about 50 cm<sup>3</sup> of 1 mol dm<sup>-3</sup> CuSO<sub>4</sub> solution.
- Place the zinc strip into a 100 cm<sup>3</sup> beaker with about 50 cm<sup>3</sup> of 1 mol dm<sup>-3</sup> ZnSO<sub>4</sub> solution.
- Use a strip of filter paper soaked in saturated potassium nitrate solution for the salt bridge
- Connect the Cu(s)|Cu<sup>2+</sup>(aq) and Zn(s)|Zn<sup>2+</sup>(aq) half-cells by connecting the metals using the crocodile clips and leads provided to the voltmeter



### Salt Bridge

The salt bridge is used to connect up the circuit. The **free moving ions** conduct the charge.

A salt bridge is usually made from a piece of filter paper (or material) soaked in a salt solution, usually **potassium nitrate**. It can also be a glass U tube containing a salt solution plugged with cotton wool

The salt should be **unreactive with the electrodes and electrode solutions**. E.g. potassium chloride would not be suitable for copper systems because chloride ions can form complexes with copper ions.


A wire is not used because the metal wire would set up its own electrode system with the solutions.

### Method

If one or both of the half cells do not contain a conducting metal, we must use an inert platinum electrode.

Set up a copper half cell using a similar arrangement to the previous one. Combine it with a Fe<sup>2+</sup>/Fe<sup>3+</sup> half-cell with a platinum electrode.

The half cell should have a mixture of acidified 1.0M iron(II) sulphate solution and an equal volume of 0.5M iron(III) sulphate solution as the electrolyte. Use a fresh salt bridge.




A platinum electrode is used because it is unreactive and can conduct electricity

Note: in the electrode system containing two solutions it is necessary to use a platinum electrode and **both** ion solutions must be of a 1M concentration so [Fe<sup>2+</sup>] = 1M and [Fe<sup>3+</sup>] = 1M .

## Core practical 13a: Redox titration

### Manganate Redox Titrations

The redox titration between  $\text{Fe}^{2+}$  with  $\text{MnO}_4^-$  (purple) is a very common exercise. This titration is self indicating because of the significant colour change from reactant to product.



The purple colour of manganate can make it difficult to see the bottom of meniscus in the burette.

If the manganate is in the burette then the end point of the titration will be the first permanent pink colour.

Colourless  $\rightarrow$  purple

#### Detailed Procedure : how much iron is in iron tablets


- Weigh accurately two 'ferrous sulfate' tablets.
- Grind up the tablets with a little 1M sulfuric acid, using a pestle and mortar.
- Through a funnel, transfer the resulting paste into a 100cm<sup>3</sup> volumetric flask. Use further small volumes of 1 M sulfuric acid to rinse the ground-up tablets into the flask.
- Then add sufficient 1M sulfuric acid to make up the solution to exactly 100cm<sup>3</sup>. Stopper the flask and shake it to make sure that all the contents are thoroughly mixed. They will not all be in solution although the  $\text{Fe}^{2+}$  ions which were present in the tablets will be dissolved.
- Titrate 10.0 cm<sup>3</sup> portions of the solution with 0.0050 M potassium manganate(VII). The end-point is marked by the first permanent purple colour.

#### Choosing correct acid for manganate titrations.

The acid is needed to supply the  $8\text{H}^+$  ions. Some acids are not suitable as they set up alternative redox reactions and hence make the titration readings inaccurate.

Only **use dilute sulfuric acid** for manganate titrations.



Insufficient volumes of sulfuric acid will mean the solution is not acidic enough and  $\text{MnO}_2$  will be produced instead of  $\text{Mn}^{2+}$ .



The brown  $\text{MnO}_2$  will mask the colour change and lead to a greater (inaccurate) volume of manganate being used in the titration.


Using a weak acid like ethanoic acid would have the same effect as it cannot supply the large amount of hydrogen ions needed ( $8\text{H}^+$ ).

It cannot be **conc HCl** as the  $\text{Cl}^-$  ions would be oxidised to  $\text{Cl}_2$  by  $\text{MnO}_4^-$  as the  $E^\circ \text{MnO}_4^-/\text{Mn}^{2+} > E^\circ \text{Cl}_2/\text{Cl}^-$



This would lead to a greater volume of manganate being used and poisonous  $\text{Cl}_2$  being produced.

It cannot be nitric acid as it is an oxidising agent. It oxidises  $\text{Fe}^{2+}$  to  $\text{Fe}^{3+}$  as  $E^\circ \text{NO}_3^-/\text{HNO}_2 > E^\circ \text{Fe}^{3+}/\text{Fe}^{2+}$




This would lead to a smaller volume of manganate being used.

#### Example 6 Manganate titration

A 2.41g nail made from an alloy containing iron is dissolved in 100cm<sup>3</sup> acid. The solution formed contains  $\text{Fe}(\text{II})$  ions.

10cm<sup>3</sup> portions of this solution are titrated with potassium manganate (VII) solution of 0.02M. 9.80cm<sup>3</sup> of  $\text{KMnO}_4$  were needed to react with the solution containing the iron.

Calculate the percentage of iron by mass in the nail.



Step1 : find moles of  $\text{KMnO}_4$

moles = conc x vol

$$0.02 \times 9.8/1000$$

$$= 1.96 \times 10^{-4} \text{ mol}$$

Step 2 : using balanced equation find moles  $\text{Fe}^{2+}$  in 10cm<sup>3</sup>

$$= \text{moles of } \text{KMnO}_4 \times 5$$

$$= 9.8 \times 10^{-4} \text{ mol}$$

Step 3 : find moles  $\text{Fe}^{2+}$  in 100cm<sup>3</sup>

$$= 9.8 \times 10^{-4} \text{ mol} \times 10$$

$$= 9.8 \times 10^{-3} \text{ mol}$$

Step 4 : find mass of Fe in  $9.8 \times 10^{-3} \text{ mol}$

$$\text{mass} = \text{moles} \times \text{RAM} = 9.8 \times 10^{-3} \times 55.8 = 0.547\text{g}$$




Step 5 : find % mass

$$\% \text{mass} = 0.547/2.41 \times 100$$




$$= 22.6\%$$

### Other useful manganate titrations

With hydrogen peroxide





#### With ethanedioate




The reaction between  $\text{MnO}_4^-$  and  $\text{C}_2\text{O}_4^{2-}$  is slow to begin with (as the reaction is between two negative ions). To do as a titration the conical flask can be heated to  $60^\circ \text{C}$  to speed up the initial reaction.

With Iron (II) ethanedioate both the  $\text{Fe}^{2+}$  and the  $\text{C}_2\text{O}_4^{2-}$  react with the  $\text{MnO}_4^-$   
 $1\text{MnO}_4^-$  reacts with  $5\text{Fe}^{2+}$  and  $2\text{MnO}_4^-$  reacts with  $5\text{C}_2\text{O}_4^{2-}$



So overall



So overall the ratio is **3**  $\text{MnO}_4^-$  to **5**  $\text{FeC}_2\text{O}_4$

#### Example 7

A 1.412 g sample of impure  $\text{FeC}_2\text{O}_4 \cdot 2\text{H}_2\text{O}$  was dissolved in an excess of dilute sulphuric acid and made up to  $250 \text{ cm}^3$  of solution.  $25.0 \text{ cm}^3$  of this solution decolourised  $23.45 \text{ cm}^3$  of a  $0.0189 \text{ mol dm}^{-3}$  solution of potassium manganate(VII).

What is the percentage by mass of  $\text{FeC}_2\text{O}_4 \cdot 2\text{H}_2\text{O}$  in the original sample?

Step 1 : find moles of  $\text{KMnO}_4$

moles = conc x vol

$$0.0189 \times 23.45/1000$$

$$= 4.43 \times 10^{-4} \text{ mol}$$

Step 2 : using balanced equation find moles  $\text{FeC}_2\text{O}_4 \cdot 2\text{H}_2\text{O}$  in  $25 \text{ cm}^3$

= moles of  $\text{KMnO}_4 \times 5/3$  (see above for ratio)

$$= 7.39 \times 10^{-4} \text{ mol}$$

Step 3 : find moles  $\text{FeC}_2\text{O}_4 \cdot 2\text{H}_2\text{O}$  in  $250 \text{ cm}^3$

$$= 7.39 \times 10^{-4} \text{ mol} \times 10$$

$$= 7.39 \times 10^{-3} \text{ mol}$$

Step 4 : find mass of  $\text{FeC}_2\text{O}_4 \cdot 2\text{H}_2\text{O}$  in  $7.39 \times 10^{-3} \text{ mol}$

$$\text{mass} = \text{moles} \times \text{Mr} = 7.39 \times 10^{-3} \times 179.8 = 1.33 \text{ g}$$

Step 5 ; find % mass


$$\% \text{mass} = 1.33/1.412 \times 100$$

$$= 94.1\%$$

## CORE PRACTICAL 13b: Redox titration

### Thiosulfate redox titration

The redox titration between  $\text{I}_2$  and thiosulfate  $\text{S}_2\text{O}_3^{2-}$  is a common exercise.



A starch indicator is added near the end point when the iodine fades a pale yellow to emphasise it.

With starch added the colour change is from blue/black to colourless


The starch should **not** be added until nearly all the iodine has reacted because the blue complex formed with high concentrations of iodine is insoluble and does not re-dissolve as more thiosulfate is added

## Core practical 14: Prepare a transition metal complex

1. Add 1.5 g of hydrated copper(II) sulfate to weighing bottle and measure the combined mass.
2. Transfer the copper (II) sulfate to a test tube and measure the mass of the empty weighing bottle.
3. Add 4 cm<sup>3</sup> of water to the test tube
4. Place the test tube in the water bath (a beaker with freshly boiled water).
5. Stir gently to dissolve the copper(II) sulfate.
6. Remove the test tube containing copper(II) sulfate solution from the water bath
7. In a fume cupboard, add 2 cm<sup>3</sup> of concentrated ammonia solution to the copper(II) sulfate solution
8. Pour the contents of the test tube into a beaker containing 6 cm<sup>3</sup> of ethanol. Stir well and cool the mixture in an ice bath.
9. Using a Buchner funnel and flask, filter the crystals. Wash the test tube with some cold ethanol and add the washings to the Buchner funnel. Finally, rinse the crystals with a little cold ethanol.
10. Scrape the crystals onto a fresh piece of filter paper and cover with a second piece of filter paper and press to dry the crystals.
11. Measure the mass of the dry crystals

The mass of the copper(II) sulfate is the difference between the two masses.

Concentrated ammonia is corrosive. Use in a fume cupboard and wear gloves



This is vacuum filtration. The apparatus is connected to a water pump which will produce a vacuum. Use if larger amounts of solid are formed.

Equation for reaction



### Loss of yield in this process

- Crystals lost when filtering or washing
- Some product stays in solution after recrystallisation
- other side reactions occurring

If the crystals are not dried properly the mass will be larger than expected which can lead to a percentage yield >100%

## Preparation of a pure organic solid and test of its purity

### Purifying an organic solid: Recrystallisation

Used for purifying aspirin


Step	Reason
1. Dissolve the impure compound in a <b>minimum volume of hot</b> (near boiling) <b>solvent</b> .	An appropriate solvent is one which will dissolve both compound and impurities when hot <b>and</b> one in which the compound itself does not dissolve well when cold. The minimum volume is used to obtain saturated solution and to enable crystallisation on cooling
2. <b>Hot filter</b> solution through (fluted) filter paper quickly.	This step will remove any insoluble impurities and heat will prevent crystals reforming during filtration
3. <b>Cool</b> the filtered solution by inserting beaker <b>in ice</b>	Crystals will reform but soluble impurities will remain in solution form because they are present in small quantities so solution is not saturated. Ice will increase the yield of crystals
4. <b>Suction filtrate</b> with a Buchner flask to separate out crystals	The water pump connected to the Buchner flask reduces the pressure and speeds up the filtration.
5 Wash the crystals with distilled water	To remove soluble impurities
6. Dry the crystals between absorbent paper	

#### Loss of yield in this process

- Crystals lost when filtering or washing
- Some product stays in solution after recrystallisation
- other side reactions occurring


If the crystals are not dried properly the mass will be larger than expected which can lead to a percentage yield >100%

Buchner flask



## 16. Core activity: The preparation of aspirin

### Detailed method for Preparation of Aspirin



Add to a 50 cm<sup>3</sup> pear-shaped flask 2.0 g of 2-hydroxybenzoic acid and 4 cm<sup>3</sup> of ethanoic anhydride.

To this mixture add 5 drops of 85% phosphoric(v) acid and swirl to mix, Fit the flask with a **reflux condenser** and heat the mixture on a boiling water bath for about 5 minutes. Without cooling the mixture, carefully add 2 **cm<sup>3</sup> of water** in one portion down the condenser.

When the vigorous reaction has ended, pour the mixture into 40 cm<sup>3</sup> of cold water in a 100 cm<sup>3</sup> beaker, stir and rub the sides of the beaker with a stirring rod necessary to induce crystallisation and, finally, allow the mixture to stand in ice bath to complete crystallisation. Collect the product by suction filtration and wash it with a little water.

### Purification stage: recrystallisation

Using a measuring cylinder, measure out 15 cm<sup>3</sup> of ethanol into a boiling tube.

Prepare a beaker half-filled with hot water from a kettle at a temperature of approximately 75 °C.

Use a spatula to add the crude aspirin to the boiling tube with ethanol and place the tube in the beaker of hot water.

Stir the contents of the boiling tube until all of the aspirin dissolves into the ethanol.

Pour the hot solution containing dissolved aspirin through a **warmed filter funnel and fluted filter** paper to hot filter

Then pour filtrate into 40 cm<sup>3</sup> of water in a conical flask.

Allow the conical flask to cool slowly and white needles of aspirin should separate.

Cool the whole mixture in an ice bath.

Filter off the purified solid under reduced pressure and allow it to dry on filter paper.

Record the mass of the dry purified solid

Aspirin is made from 2-hydroxybenzoic acid which contains a phenol group. In the reaction the phenol group is turned into an ester by reacting it with the reactive ethanoic anhydride

Ethanoic anhydride is used instead of acid chlorides because it is cheaper, less corrosive, less vulnerable to hydrolysis, and less dangerous to use.

The excess ethanoic anhydride will hydrolyse and the contents of the flask will boil.

Avoid naked flames due to flammability of ethanol

This step will remove any insoluble impurities and heat will prevent crystals reforming during filtration

Soluble impurities will remain in solution form because they are present in small quantities so solution is not saturated. Ice will increase the yield of crystals

### Detailed method for Nitration Procedure

Measure 2.5 cm<sup>3</sup> of methyl benzoate into a small conical flask and then dissolve it in 5 cm<sup>3</sup> of concentrated sulfuric acid. When the liquid has dissolved, cool the mixture in ice.

Prepare the nitrating mixture by adding **drop by drop** 2 cm<sup>3</sup> of concentrated sulfuric acid to 2 cm<sup>3</sup> of concentrated nitric acid. **Cool this mixture in ice** as well.

Now add the nitrating mixture drop by drop from a dropping pipette to the solution of methyl benzoate. Stir the mixture with a thermometer and keep the temperature **below 10 °C**. When the addition is complete, allow the mixture to stand at room temperature for another 15 minutes.

After this time, pour the reaction mixture on to about 25 g of crushed ice and stir until all the ice has melted and crystalline methyl 3-nitrobenzoate has formed.

Then use same purification method as in aspirin above

Conc acids are corrosive- wear gloves

The acids react together to make the NO<sub>2</sub><sup>+</sup> ion

This reaction is exothermic so acids are kept cool and acid is added dropwise

The temperature is kept low at this stage to prevent multiple substitution of nitro groups on the benzene ring

## Measuring melting point

One way of testing for the degree of purity is to determine the melting "point", or melting range, of the sample.


If the sample is very pure then the melting point will be a sharp one, at the same value as quoted in data books.

If **impurities** are present (and this can include solvent from the recrystallisation process) the **melting point will be lowered** and the sample will **melt over a range** of several degrees Celsius

Melting point can be measured in an electronic melting point machine or by using a practical set up where the capillary tube is strapped to a thermometer immersed in some heating oil.

In both cases a small amount of the sample is put into a capillary tube. The tube is heated up and is **heated slowly near the melting point**

Comparing an experimentally determined melting point value with one quoted in a data source will verify the degree of purity.




Sometimes an error may occur if the temperature on the thermometer is not the same as the temperature in the actual sample tube.

## Steam distillation

In steam distillation steam is passed into the mixture and the product vapour is distilled off with the water and condensed

### Advantage of steam distillation:

The product distils at a lower temperature which can prevent decomposition of the product if it has a high boiling point



## Solvent extraction

Mix organic solvent and oil-water mixture in a separating funnel then separate the oil layer.  
Distil to separate oil from organic solvent  
Add anhydrous  $\text{CaCl}_2$  to clove oil to dry oil  
Decant to remove  $\text{CaCl}_2$

Separating funnel



## Separation of species by thin-layer chromatography

### Method: Thin-layer chromatography

- Wearing gloves**, draw a **pencil line** 1 cm above the bottom of a TLC plate and mark spots for each sample, equally spaced along line.
- Use a capillary tube to add a **tiny drop** of each solution to a different spot and allow the plate to air dry.
- Add solvent to a chamber or large beaker with a lid so that is no more than **1cm in depth**
- Place the TLC plate into the chamber, **making sure that the level of the solvent is below the pencil line**. Replace the lid to get a **tight seal**.
- When the level of the solvent **reaches about 1 cm from the top of the plate**, remove the plate and mark the solvent level with a pencil. Allow the plate to **dry in the fume cupboard**.
- Place the plate under a **UV lamp** in order to see the spots. Draw around them lightly in pencil.
- Calculate the  $R_f$  values of the observed spots.

Wear plastic gloves to prevent contamination from the hands to the plate

**pencil line** –will not dissolve in the solvent

**tiny drop** – too big a drop will cause different spots to merge

**Depth** of solvent– if the solvent is too deep it will dissolve the sample spots from the plate


**lid**– to prevent evaporation of toxic solvent

Will get more accurate results if the solvent is allowed to rise to near the top of the plate but the  $R_f$  value can be calculated if the solvent front does not reach the top of the plate

dry in a **fume cupboard** as the solvent is toxic

**UV lamp** used if the spots are colourless and not visible

$$R_f \text{ value} = \frac{\text{distance moved by amino acid}}{\text{distance moved by the solvent}}$$



Separation by chromatography depends on the balance between solubility in the moving phase and retention in the stationary phase.

A solid stationary phase separates by adsorption,  
A liquid stationary phase separates by relative solubility


If the stationary phase was polar and the moving phase was non-polar e.g. Hexane. Then non-polar compounds would pass through the plate more quickly than polar compounds as they would have a greater solubility in the non-polar moving phase. (Think about intermolecular forces)

$R_f$  values are used to identify different substances.



## Practical: Working out equilibrium constant Kc

A common experiment is working out the equilibrium constant for an esterification reaction. Ethanol and ethanoic acid are mixed together with a sulfuric acid catalyst.



### Method

#### Part 1 Preparing the equilibrium mixture

- 1 Use burettes to prepare a mixture in boiling tube of carboxylic acid, alcohol, and dilute sulfuric acid.
- 2 Swirl and bung tube. Leave the mixture to reach equilibrium for one week

#### Part 2 Titrating the equilibrium mixture

- 1 Rinse a 250 cm<sup>3</sup> volumetric flask with distilled water. Use a funnel to transfer the contents of the boiling tube into the flask. Rinse the boiling tube with water and add the washings to the volumetric flask.
- 2 Use distilled water to make up the solution in the volumetric flask to exactly 250 cm<sup>3</sup>. Stopper the flask, then invert and shake the contents thoroughly.
- 3 Use the pipette to transfer 25.0 cm<sup>3</sup> of the diluted equilibrium mixture to a 250 cm<sup>3</sup> conical flask.
- 4 Add 3 or 4 drops of phenolphthalein indicator to the conical flask.
- 5 Set up the burette with sodium hydroxide solution..
- 6 Add the sodium hydroxide solution from the burette until the mixture in the conical flask just turns pink. Record this burette reading in your table.
- 7 Repeat the titration until you obtain a minimum of two concordant titres.

The sodium hydroxide will react with the sulfuric acid catalyst and any unreacted carboxylic acid in the equilibrium mixture

There are many different calculations that can be based on this experiment. Let's look at general stages. Not all calculations will use all the stages.

#### Working out initial amount of moles of reactants

The amount of moles of alcohol and carboxylic acid can be calculated from the densities and volumes of liquids added

Mass = density x volume

then

Moles = mass / Mr

The initial amount of moles of acid catalyst used is usually determined by titrating a separate sample of catalyst with sodium hydroxide

#### Working out equilibrium amount of moles of acid present from the titre results

39.0 cm<sup>3</sup> of 0.400 mol dm<sup>-3</sup> sodium hydroxide was used in the above titration. The initial moles of sulphuric acid was 5x10<sup>-4</sup> mol. Calculate the moles of ethanoic acid present at equilibrium

Amount of NaOH = vol X conc  
= 0.039 x 0.400  
= 0.0156 mol

So total amount of H<sup>+</sup> present in 25cm<sup>3</sup> = 0.0156 mol

So total amount of H<sup>+</sup> present in 250cm<sup>3</sup> = 0.156 mol

Total mol acid present = moles of carboxylic acid + moles of acid catalyst

So

Amount of carboxylic acid at equilibrium = 0.156 - (5x10<sup>-4</sup> x 2)  
= 0.155 mol

X 2 because H<sub>2</sub>SO<sub>4</sub>  
has 2 H<sup>+</sup>

The pink colour of the phenolphthalein in the titration can fade after the end-point of the titration has been reached because the addition of sodium hydroxide may make the equilibrium shift towards the reactants

### Working out equilibrium amount of moles of other substances

Calculate the equilibrium amount of ethanol, ethyl ethanoate and water if there were initially 0.400 mol of ethanol and 0.500 mol of ethanoic acid and at equilibrium there were 0.155 mol of ethanoic acid.

$$\begin{aligned}\text{Amount of ethanoic acid that reacted} &= \text{initial amount} - \text{equilibrium amount} \\ &= 0.5 - 0.155 \\ &= 0.344 \text{ mol}\end{aligned}$$

$$\begin{aligned}\text{Amount of ethanol at equilibrium} &= \text{initial amount} - \text{amount that reacted} \\ &= 0.400 - 0.344 \\ &= 0.056 \text{ mol}\end{aligned}$$

$$\begin{aligned}\text{Amount of ethyl ethanoate at equilibrium} &= \text{initial amount} + \text{amount that formed} \\ &= 0 + 0.344 \\ &= 0.344 \text{ mol}\end{aligned}$$

$$\begin{aligned}\text{Amount of water at equilibrium} &= \text{initial amount} + \text{amount that formed} \\ &= 0 + 0.344 \\ &= 0.344 \text{ mol}\end{aligned}$$

The amount of water at equilibrium would not really be 0 as there would be water present in the acid catalyst

### Calculating the equilibrium constant

Finally calculate the equilibrium constant.

To work out equilibrium concentrations divide the equilibrium amounts by the total volume. Then put in Kc expression

$$K_c = \frac{[\text{CH}_3\text{CO}_2\text{CH}_2\text{CH}_3] [\text{H}_2\text{O}]}{[\text{CH}_3\text{CO}_2\text{H}] [\text{CH}_3\text{CH}_2\text{OH}]}$$

In order to confirm that one week was sufficient time for equilibrium to be established in the mixture from Part 1, several mixtures could be made and left for different amount of time. If the resulting Kc is the same value then it can be concluded the time is sufficient

## Spectrophotometry

If visible light of increasing frequency is passed through a sample of a coloured complex ion, some of the light is absorbed.

The amount of light absorbed is proportional to the concentration of the absorbing species (and to the distance travelled through the solution).

Some complexes have only pale colours and do not absorb light strongly. In these cases a suitable ligand is added to intensify the colour.

Absorption of visible light is used in spectrometry to determine the concentration of coloured ions.


### method

- Add an appropriate ligand to **intensify** colour
- Make up solutions of known concentration
- Measure absorption or transmission
- Plot graph of absorption vs concentration
- Measure absorption of unknown and compare

Spectrometers contain a coloured filter. The colour of the filter is chosen to only allow the wavelengths of light through that would be most strongly absorbed by the coloured solution.

### Detailed method- measuring absorption of copper solutions

- Take nine 100cm<sup>3</sup> graduated flasks and pipette 20cm<sup>3</sup> of 2 mol dm<sup>-3</sup> ammonia solution into each one. Use the 0.05M solution of aqueous copper sulfate to make up solutions which are concentrations of 0.005 to 0.04 mol dm<sup>-3</sup> of [Cu(NH<sub>3</sub>)<sub>4</sub>(H<sub>2</sub>O)<sub>2</sub>]<sup>2+</sup>
- Mix each solution thoroughly.
- Insert the red filter into the colorimeter.
- Use a cuvette with distilled water to zero the colorimeter
- Then put each prepared solution in cuvette and measure the absorbance of each solution.
- Plot graph of absorption vs concentration
- Measure absorption of unknown solution and determine its concentration from the calibration curve




0.005 to 0.04 mol dm<sup>-3</sup> solutions of [Cu(NH<sub>3</sub>)<sub>4</sub>(H<sub>2</sub>O)<sub>2</sub>]<sup>2+</sup>




colorimeter

Calibration curve



## Autocatalytic Reaction between Ethanedioate and Manganate ions

Potassium manganate(VII) will oxidise ethanedioic acid (oxalic acid) to carbon dioxide and water, in the presence of an excess of acid:




### Detailed method

1. Prepare a reaction mixture according to the table, using measuring cylinders.
2. Some members of your group should use mixture 1 and some mixture 2. The results should then be shared.
3. Add 50 cm<sup>3</sup> of 0.02 M potassium manganate(VII) and start timing. Shake the mixture for about half a minute to mix it well.
4. After about a minute use a pipette to withdraw a 10.0 cm<sup>3</sup> portion of the reaction mixture and empty into a conical flask.
5. Note the time and add about 10cm<sup>3</sup> of 0.1 M potassium iodide solution. This stops the reaction and releases iodine equivalent to the remaining manganate(VII) ions.
6. Titrate the liberated iodine with 0.01 M sodium thiosulfate, adding a little starch solution near the end-point. Record the titre of sodium thiosulfate.
7. Remove further portions every 3 or 4 minutes and titrate them in the same way. Continue until the titre is less than 3 cm<sup>3</sup>.



Solution	Mixture 1	Mixture 2
0.2 M ethanedioic acid	100cm <sup>3</sup>	100cm <sup>3</sup>
0.2 M manganese(II) sulfate	--	15cm <sup>3</sup>
1 M sulfuric acid	10 cm <sup>3</sup>	10cm <sup>3</sup>
Water	90 cm <sup>3</sup>	75 cm <sup>3</sup>

### Explanation of results

The autocatalysis by Mn<sup>2+</sup> in titrations of C<sub>2</sub>O<sub>4</sub><sup>2-</sup> with MnO<sub>4</sub><sup>-</sup>



#### Catalysed alternative route




The initial uncatalysed reaction is **slow** because the reaction is a collision between **two negative ions** which **repel each other** leading to a **high activation energy**.

The Mn<sup>2+</sup> ions produced act as an **autocatalyst** and therefore the reaction starts to speed up because they bring about the alternative reaction route with lower activation energy.

The reaction eventually slows as the **MnO<sub>4</sub><sup>-</sup> concentration drops**.

This is an example of **autocatalysis** where one of the products of the reaction can catalyse the reaction.



### Alternative method for following the reaction rate

This experiment can be done by removing samples at set times and titrating to work out the concentration of MnO<sub>4</sub><sup>-</sup>. It could also be done by use of a spectrometer measuring the intensity of the purple colour. This method has the advantage that it **does not disrupt the reaction mixture**, using up the reactants and it leads to a much **quicker determination of concentration**.